# **General Info**



## Finally a walking cane that helps you walk better! Specifically engineered to meet the needs of everyday users.

www.fizio.ca

Increases stability Improves walking Improves your posture Reduces pain Ultra-lightweight Stands on its own! Unique design Improves self image! Available in 6 colours Silver - Black - Bronze - Blue - Burgundy - Red

*Covered by most insurances* 

Adjustable height: Standard model: 29.5"to 35.5" Short model: 26" to 29" 250 Lbs Maximum User weight capacity



Rubber Sole Series 100



Rubber Sole Left -Series 300



Rubber Sole Right -Series 300

### To order, Call us toll free at 1-877-424-4114 We deliver everywhere in Canada



**Evolution in motion** 

info@FIZIO.ca / Tel : (450) 424-9869 / Fax : (450) 510-0842

## **Rubber soles**

# Fizio<sup>™</sup> Walking Stabilizer

## **Finally a walking cane that helps you walk better!** Specifically engineered to meet the needs of everyday users.

# Rubber soles Fully interchangeable

Series 100



**Dimensions :** 2.25"wide 4.5"long 1.5"high **Black** 60g

#### **Specifications**

- Rocker bottom. Improves gait pattern.
- Good for all seasons. No need of an ice pick.
- Rubber lips create a suction effect on winter and slippery surfaces.
- Grooves to improve grip on snowy surfaces.
- Rubber lips makes the FIZIO stand on its own.
- Made of anti-slip rubber compound.

#### Who is it for?

- For people who need a slight walking support.
- For people who are still mobile, who can walk for several minutes during the day.

Series 300





Dimensions: 4"wide x 4.5"long

x 2.25"high / Black / 230g

Specifications - More lateral stability.

- Can replace a quad cane.
- Fits on stairs.
- Helps people re-learn how to walk.
- Flat surface at bottom.
- Rocker bottom removed.
- Made of anti-slip rubber compound.

#### Who is it for?

- For people who are unstable and want more security.
- For people who put a lot of weight on their cane.
- For people who already use a quad cane.
- For people who had a surgery or stroke.
- For people in rehabilitation.

#### To order, Call us toll free at 1-877-424-4114 We deliver everywhere in Canada

**Evolution in motion** 

www.fizio.ca

# Fizio<sup>™</sup> Walking Stabilizer

# **Additional Info**

### **Finally a walking cane that helps you walk better!** Specifically engineered to meet the needs of everyday users.



#### Use it the handle FACING FORWARD



# Coverage by most insurances

Yes, most private and corporate insurance companies will pay up to 80% of the cost of the FIZIO walking stabilizer. You will need to provide them with a prescription from a specialist such as physiatrist, orthopedic surgeon, neurologist, rheumatologist or geriatrician and the invoice from the purchase of your FIZIO<sup>™</sup> walking stabilizer. It is important that the doctor specifically writes on the prescription **"FIZIO<sup>™</sup> walking stabilizer".**  How to use the ROCKER BOTTOM:

Put it down heel first...



...then let it roll ...




...and push off with the toe.



You will find that, unlike a cane, you will have normal arm movement as you walk, increasing your comfort and causing less fatigue in your arm and wrist.

#### What MODEL is right for you?

SHORT : 26" TO 29" STANDARD : 29.5" TO 35.5"



Stand tall, with your shoes on, your arms straight falling alongside of your body. The recommended height is the distance from the floor to the beginning of your wrist, at the line where your wrist bends.

### To order, Call us toll free at 1-877-424-4114 We deliver everywhere in Canada



**Evolution in motion** info@FIZIO.ca / Tel : (450) 424-9869 / Fax : (450) 510-0842

# **Price List**

# Fizio<sup>™</sup> Walking Stabilizer

**Finally a walking cane that helps you walk better!** Specifically engineered to meet the needs of everyday users.



# www.fizio.ca

To order, Call us toll free at 1-877-424-4114 Free delivery everywhere in Canada



## **Order Form**

# Fizio<sup>™</sup> Walking Stabilizer

Please check the boxes below and sign the form to complete your order

Packages:	(Free Delivery)	Price:		color:	Side:	Height:
Simple Package 1 FIZIO, Rubber sole series 100 included		\$285	Bronze	Blue Burgundy		Short Standard
Complete Package     1 FIZIO, Rubber sole series 100 included     1 Rubber sole series 100, 1 Wrist strap			Silver	Blue Burgundy Red	1	Short
Simple Stability Package 1 FIZIO, Rubber sole series 300 included		\$305	Silver	Blue Burgundy Red	└ □ Right □ Left	Short
Complete Stabili     1 FIZIO, Rubber sole     1 Rubber sole series	\$315	Silver	Blue Burgundy Red	□ Right □ Left I	Short	
	harges extra ht with a package)	ing charges	6			
Wrist strap \$9         Rubber sole, series 100 \$20         Rubber Sole Series 200 Right \$35         Rubber Sole Series 200 Left \$35						
				VWW	<i>ı</i> .fiz	cio.ca
Billing Address						
LAST NAME			Γ NAME			
ADDRESS						
CITY			PROVINCE POSTAL CODE			
TELEPHONE #		FAX #	OR EMAIL			
Shipping Address	Same as Billing address					
ADDRESS						
CITY		PRO\	INCE		POSTAL CODE	
Payment information			Credit Card:		54	Master Card
Sub-Total: \$			Card Holder			
Shipping charges: \$			name: _			
			Card Number:			
Total amount to be put on my credit card: \$		_	Expiration	/		
I agree to pay the total amount above		e :	Date: _			
						ORM TO:
Signature:			Date	450	)-510-	-0842

#### MEDICAL RECOMMENDATION

NAME OF INSTITUTION	NATURE OF RECOMMENDATION
PROGRAM	
ADDRESS	
CITY	
PROVINCE POSTAL CODE	

#### (PATIENT INFORMATION )

LAST NAME	FIRST NAME
ADDRESS	
CITY	PROVINCE POSTAL CODE
TELEPHONE #	FAX #

FILE #	WSIB	ACCIDENT DATE:	YEAR	MONTH	DAY
NAME OF INDEMNISATION AGENT					
TELEPHONE #	FAX #				

(PATIENT'S DISABILITY)				
RECOMMENDED EQUIPMENT				
Advise if equipment is to be worn daily			Recommended time frame	
U YES			MONTHS	
(MEDICAL IMPROVEMENT SOUGHT )				
	GREATER AUT	ONOMY	IMPROVE POSTURE	
REDUCE RISK OF FALLS	REDUCE PAIN		IMPROVE QUALITY OF LIF	E
		UE		

RECOMMENDED BY	TITLE	TELEP	HONE #			
SIGNATURE			DATE	YEAR	MONTH	DAY